

# The Gripper System











## LESS HOLDING, MORE SURGERY

Human capital is vital in the OR. The Gripper™ System ensures stable and flexible retractor positioning, allowing nursing staff to be effectively deployed on more valuable tasks.

## 19 02 v3 ME GE Gripper - © MedEnvision BV 2024 - All rights reserv

### Gripper™ System

THE ASSISTANT IN A BOX

The Gripper™ System allows the surgeon to control the retractor without an assistant

A well-managed OR can boost surgical turnover, reduce complications, enhance patient outcomes and satisfaction, improve ergonomics, and increase productivity.

The Gripper™ System puts hip, knee and shoulder surgeons back in charge of their instruments. They determine the positioning, the force with which the retractors are placed and how long they stay in place. This allows the surgeons to be 100% focused on their job while the scrub is always one step ahead because he/she doesn't have to hold retractors.



#### 24/7 AVAILABILITY

The Gripper™ System is the assistant that is always sterile and available 24/7. In case of limited availability of staff, the procedure will not have to be postponed and will therefore not affect the efficient functioning of the OR.



#### **EFFICIENT STAFF ALLOCATION**

Because the Gripper™ System replaces hand held retractors, they enable nurses to be allocated to more productive work and urgent operations. This enables parallel processing and efficient allocation of the OR staff.



#### BETTER TISSUE PRESERVATION

By using the Gripper<sup>TM</sup> System the surgeon determines the force applied on the retractor which reduces the cutting effect on the bone and soft tissues.

During fluoroscopy the retractor does not need to be removed which again reduces damage to the soft tissues.



#### **INCREASED STABILITY**

The surgeon places the Gripper™ System with the correct force and direction so that the retractor can be positioned accurately and stable. The Gripper™ System remains stable in its position even under the most demanding circumstances.



#### **OPTIMAL VISIBILITY**

The Gripper<sup>TM</sup> System is connected to the side of the table offering a project view on the surgical field and anatomical landmarks, not obstructed by staff holding retractors.



#### **IMPROVED EFFICIENCY**

By using the Gripper™ System the surgeon can focus on his surgery. The Gripper™ System gives the surgeon more control over the procedure and contributes to standardizing the process. Standardization leads to optimized efficiency and increased productivity.

Item	Description	Art. code	Unit	Packed per
Gripper™	Self Retracting Device	ME-GR1000	PC	10
Fixation Devices	Fixation Device Hip Set (regular table)	ME-FD1000	SET	1
	Fixation Device Table Adapter (traction table)	ME-FD3000	SET	1
	Fixation Device Knee Set	ME-FD4000	SET	1
	Fixation Device Shoulder Set	ME-FD5000	SET	1

